


Charge transfer as a diagnostic of excitation dynamics in a MOT

G. Veshapidze, M. L. Trachy, H. A. Camp,
M. A. Gearba, M. H. Shah, X. Flechard, H. U. Jang,
H. Nguyen, J. Blicck, R. Bredy, C. W. Fehrenbach,
B. D. DePaola


J.R. Macdonald Laboratory
Department of Physics, Kansas State University
Manhattan, Kansas

Charge Transfer with MOTRIMS

$$Q = E_i - E_f$$
$$= -V_p \boxed{P_{||}} - \frac{1}{2} m_e V_p^2$$




Apparatus




Ref: R. Bredy *et al.*, Nucl. Instrum. and Meth. Phys. Res. B 205, 191 (2003)

Charge Transfer with MOTRIMS



$$Q = E_i - E_f$$




$$A_i = cn_i\sigma_i$$

A_i = area under the peak

n_i = number of atoms in the i -th state

σ_i = Charge transfer cross section

Diagnostics with Charge Transfer



$$f = \frac{n_p}{n_p + n_s}$$

$$= \frac{A_p}{A_p + A_s \frac{\sigma_p}{\sigma_s}}$$


$$\frac{\sigma_p}{\sigma_s} = 11.29 \pm 0.66$$

$$A_s = c n_s \sigma_s$$

$$A_p = c n_p \sigma_p$$

Ref: X. Flechard *et al.*, PRL 91, 243005 (2003)

Diagnostics with Charge Transfer



Trapping Laser is off
only for a short time, so that

$$n_{tot} = \text{Const}$$

$$\Delta n_p = -\Delta n_s$$


$$\Delta A_p / \sigma_p = -\Delta A_s / \sigma_s$$

$$\sigma_p / \sigma_s = -\Delta A_p / \Delta A_s$$

Ref: X. Flechard *et al.*, PRL 91, 243005 (2003)

Application Example: Excited State Fraction

$$f \equiv \frac{n_{ext}}{n_{tot}}$$



- *Two level system*
- *A single plain wave*
- *Linearly polarized*
- *Intensity is low*

$$f = \frac{I/I_s}{1 + 2I/I_s + (2\delta/\Gamma)^2}$$

$$I_s = \frac{2\pi hc\Gamma}{3\lambda^3}$$



Ref: W. Demtröder, Laser Spectroscopy. (Springer, 2002)

Application Example: Excited State Fraction

$$f \equiv \frac{n_{ext}}{n_{tot}}$$

- Absolute Photo-Ionization Cross Sections
- Cold Atoms Collisions Cross Sections
- Photo-Association Rates
- Total Number of Atoms in a MOT

The Problem




- *Multilevel system*
- *Six plane waves*
- *Circularly polarized*
- *Intensity is not low*

$$f \stackrel{?}{=} \frac{I/I_s}{1 + 2I/I_s + (2\delta/\Gamma)^2}$$

$$I_s \stackrel{?}{=} \frac{2\pi h c \Gamma}{3\lambda^3}$$

Results: Simple Model (1-Parameter) Fit



$$f = \frac{I/I_s}{1 + 2I/I_s + (2\delta/\Gamma)^2}$$

$$s = \frac{I/I_s}{1 + (2\delta/\Gamma)^2}$$

$$f = \frac{s}{1 + 2s}$$

Theory says:


$$I_s = 3.28 \text{ mW/cm}^2$$

1-Parameter fit
gives:


Ref: M. Shah *et al.*, Phys. Rev. A 75, 053418 (2007) $I_s = 9.2 \text{ mW/cm}^2$

Application example: STIRAP


- STIRAP = Stimulated Raman Adiabatic Passage



Efficiently move the population from g to e without passing through i



= Intuitive Pulse Ordering




= Counter-Intuitive Pulse Ordering

Applications of STIRAP


- State Preparation in Quantum Computing
- Reaction Channel Activation in Chemistry
- Atom Optics
- High-Order Harmonics Generation

Our System for STIRAP Studies: ^{87}Rb



Ref: Gearba *et al.*, PRA 76, 013406 (2007)

Results




Counter-Intuitive

Intuitive

Ref: Gearba *et al.*, PRA 76, 013406 (2007)

Acknowledgements

Team MOTRIMS



Hyounguk Jang, Mudessar Shah, Brett DePaola, Marc Trachy, Giorgi Veshapidze

Special Thanks to
MOTRIMS Team
JRM Labs People
Thanks for your attention!